

## Clinical Evaluation Results Prove That Bioenteric Is More Sensitive Than Standard Culture

### Waterford Regional Hospital, HSE South, Ireland

Pathogen	Standard Culture	BioEnteric
Salmonella	4	5*
E coli 0157	0	3*
Shigella	0	0
Campylobacter	17	25*
Negatives	814	802
Samples Tested	835	835

### Mid Western Regional Hospital, HSE West, Ireland

Pathogen	Standard Culture	BioEnteric
Salmonella	3	4*
E coli 0157	2	5*
Shigella	0	1*
Campylobacter	19	25*
Negatives	526	515
Samples Tested	550	550

\*Discrepancies confirmed using alternative PCR method